OmniCeram évolution flow 2.0

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Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Etching Gel and subsequent application of a light cure bonding agent (e.g. OmniCite 2.0) is recommended. Otherwise use a self etching bonding agent (e.g. OmniCite All'N'One) before application of OmniCeram évolution flow 2.0.

Placement of OmniCeram évolution flow 2.0

Apply OmniCeram évolution flow 2.0 in the selected shade. Place it by using a suitable instrument. Transparent matrix strips may be used.

Application by syringe:

Screw off the cap of the syringe and substitute it by e.g. OmniCeram évolution flow Tips. After use discard the Tips and attach the cap to the syringe again. Tips are for single use only due to hygienic reasons.

Curing of OmniCeram évolution flow 2.0

For an optimum result apply OmniCeram évolution flow 2.0 in **layers of max. 1 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity of at least 1000mW/cm²** and a wavelength range of 400 to 500 nm (blue light) cure each increment as follows:

Lighter shades (e.g. A1, A2, A3, B1, Bleach) **20 seconds**

Darker shades (e.g. A3.5) **30 seconds**

Hold the light emission window as close as possible to the filling material

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings

- PL-Tips are for single use only due to hygienic reasons.
- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of the composite restorative.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.
- Keep away from children!

Composition

Dimethacrylates, dental glass, silicon dioxide, photo initiators.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

SanaPro Dental GmbH warrants this product will be free from defects in material and manufacture. SanaPro Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and SanaPro Dental GmbH's sole obligation shall be repair or replacement of the SanaPro Dental GmbH product.

Limitation of Liability

Except where prohibited by law, SanaPro Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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