OMNI DENT®

PRAXIS

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We care. You smile.

Content: (→) 50 ml (70g) Etching Gel,
10 x Applicator syringes,
2 x Adapter, 50 x Applicator Tips.

Inhoud: (→) 50 ml (70g) Etsgel spuit,
10 x Applicatiecanules,
2 x Adapter, 50 x Applicatie-Tips.

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CE 0124

Safety guidelines:

H 314 Causes severe skin burns and eye damage. Wear protective gloves/eye protection. IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a long time with water. Remove any contact lenses if possible. Continue to rinse. Call a POISON CENTRE or a doctor/physician immediately. IF SWALLOWED: Rinse the mouth. DO NOT induce vomiting. IN THE EVENT OF CONTACT WITH THE SKIN: Wash the area with a lot of soap and water. Store in a closed container. Dispose of the contents/container in accordance with the local/regional/national/international provisions. Contains orthophosphoric acid.

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Product description:

Etching gel is a 37% phosphoric acid gel with optimum thixotropic flow properties. It is used for etching dental substances during enamel etching procedures and total-etch techniques. If the safety notes and instructions on use are observed, the product has no side-effects according to the current state of knowledge. Underfilling materials can be applied by the etching gel with a longer application time.

- Only for dental use!
- Do not use after the use-by date has expired.
- Close the syringe immediately after use.

Instructions for use:

Filling empty application syringes:

An empty application syringe is pulled towards the supply syringe using the adaptor (turn 180° clockwise). The small syringe is filled by pressing the large syringe plunger down carefully and slowly.

The small syringe plunger may also be pulled up slightly at the same time, and you may have to make sure that the small syringe is not overfilled or the plunger pressed out of the syringe.

Then the small application syringe is twisted off the adaptor

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e della soml syringe is closed tightly again. A cannula is placed onto the application syringe.

Preparation of filled application syringes:

A filled application syringe is taken out of the packaging and the cap is removed. An application cannula is screwed onto the syringe in a clockwise direction and tightened. The etching gel can now be easily applied by carefully pressing the syringe plunger.

Preparation:

The area of the dental substance to be etched must be clearly marked.

- 1.) The use of a rubber dam is recommended! Any contact with soft tissue (skin, oral mucosa, eyes etc.) is to be avoided. If etching gel gets into the soft tissue, rinse the area immediately with plenty of water and, if necessary, contact a doctor.
- 2.) Before preparation, clean the tooth using a rubber cup with fluoride and fat-free polishing paste.
- 3.) Remove cleaning residue carefully with water and clean and dry the surface of the tooth with oil-free compressed air.
- 4.) Prepare the cavity accordingly, bevel the melted edges to achieve permanent enamel junctions. Rinse the cavity thoroughly with water and dry with oil-free compressed air. Avoid overdrying!

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5.) Protect dental pulp in deep cavities. (e.g. cover with a thin layer of suitable underfilling material). The cavity should be dry and free of contamination, whereby the overdrying of the dentine is to be prevented.

Etching process

6.) Apply the etching gel to the melted edges of the cavity and, where necessary, to the dentine.

Recommended application time:

- 15 - 30 sec. on dental enamel
- 10 - 15 sec. on dentine
- 50 - 80 sec. on milk teeth and
- Teeth with high fluoride content
- 7.) Rinse for at least 20 sec. with water.
- 8.) Dry with oil and water-free compressed air. The etched surface is not allowed to be touched until the treatment is started, or the surface has been contaminated by saliva. However, should contamination still occur after the etching and drying, the area must be etched again briefly (approx. 10 sec.).

The etched enamel should look chalky white!

After the etching, you can continue with the application of adhesives etc. in accordance with the manufacturer's instructions for use.

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